

Alexa Fluor 647 anti-mouse IFN- γ

Catalog No: #FC1309

Package Size: #FC1309-1 25Tests #FC1309-2 100Tests

Description

Product Name	Alexa Fluor 647 anti-mouse IFN- γ
Host Species	rat
Clonality	Monoclonal
Clone No.	s-XMG1.2
Isotype	rat IgG1
Applications	FC
Species Reactivity	mouse
Conjugates	AF Dye
Target Name	IFN- γ
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.