

# Purified anti-human CD4

Catalog No: #FC1427

Package Size: #FC1427-1 100ug #FC1427-2 1mg

## Description

|                    |                         |
|--------------------|-------------------------|
| Product Name       | Purified anti-human CD4 |
| Host Species       | mouse                   |
| Clonality          | monoclonal              |
| Clone No.          | s-OKT4                  |
| Species Reactivity | human                   |
| Conjugates         | Unconjugated            |
| Storage            | Store at 4°C.           |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.