

Purified anti-human CD45

Catalog No: #FC1488

Package Size: #FC1488-1 100ug #FC1488-2 1mg



Signalway Antibody

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                    |                          |
|--------------------|--------------------------|
| Product Name       | Purified anti-human CD45 |
| Host Species       | mouse                    |
| Clonality          | monoclonal               |
| Clone No.          | s-2D1                    |
| Species Reactivity | human                    |
| Conjugates         | Unconjugated             |
| Storage            | Store at 4°C.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.