## Purified anti-human HLA-AoO BoO CoO HLA Class IoO

Catalog No: #FC1550

Package Size: #FC1550 100 ug



Support: tech@signalwayantibody.com

## Description

Product Name	Purified anti-human HLA-AoO BoO CoO HLA Class IoO
Host Species	Mouse
Clonality	monoclonal
Clone No.	FM1001
Applications	IHC,IF,FC
Species Reactivity	human
Conjugates	Unconjugated
Storage	Store at 4°C.

## **Application Details**


IHC: 1:50-1:200

IF: 1:50-1:200

FC:5 ΞΌI/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 106 cells can be tested per experiment. Refer to the detailed scenario as you test.

## **Images**



Formalin-fixed, paraffin-embedded Human breast carcinoma tissue stained for human HLA-AoO BoO C using FC1550 at 1/50 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.