

Purified anti-human HLA-DP (MHC Class II)
Monomorphic

Catalog No: #FC1552

Package Size: #FC1552-1 100ug #FC1552-2 1mg

Description

Product Name	Purified anti-human HLA-DP (MHC Class II) Monomorphic
Host Species	mouse
Clonality	monoclonal
Clone No.	s-B7/21
Species Reactivity	human
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.