Product Datasheet

Purified anti-human/mouse IL-6

Catalog No: #FC1563

Package Size: #FC1563-1 100ug #FC1563-2 1mg



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human/mouse IL-6
Host Species	human
Clonality	monoclonal
Clone No.	s-Hu137
Species Reactivity	human/mouse
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 106 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.