

Purified anti-Mouse/Human CD49d

Catalog No: #FC1599

Package Size: #FC1599-1 100ug #FC1599-2 1mg

Description

Product Name	Purified anti-Mouse/Human CD49d
Host Species	rat
Clonality	monoclonal
Clone No.	s-PS/2
Species Reactivity	Mouse/Human
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.