

# Purified anti-mouse CD70

Catalog No: #FC1604

Package Size: #FC1604-1 100ug #FC1604-2 1mg

## Description

|                    |                          |
|--------------------|--------------------------|
| Product Name       | Purified anti-mouse CD70 |
| Host Species       | mouse                    |
| Clonality          | monoclonal               |
| Clone No.          | s-TAN1-7                 |
| Species Reactivity | mouse                    |
| Conjugates         | Unconjugated             |
| Storage            | Store at 4°C.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.