

Purified anti-mouse CD90

Catalog No: #FC1610

Package Size: #FC1610-1 100ug #FC1610-2 1mg



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse CD90
Host Species	rat
Clonality	monoclonal
Clone No.	s-HK2.1
Species Reactivity	mouse
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.