

## Purified anti-mouse CD178

Catalog No: #FC1631

Package Size: #FC1631-1 100ug #FC1631-2 1mg

## Description

|                    |                           |
|--------------------|---------------------------|
| Product Name       | Purified anti-mouse CD178 |
| Host Species       | Armenian Hamster          |
| Clonality          | monoclonal                |
| Clone No.          | s-MFL3                    |
| Species Reactivity | mouse                     |
| Conjugates         | Unconjugated              |
| Storage            | Store at 4°C.             |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.