

Purified anti-mouse CXCL9

Catalog No: #FC1653

Package Size: #FC1653-1 100ug #FC1653-2 1mg



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse CXCL9
Host Species	Armenian Hamster
Clonality	monoclonal
Clone No.	s-MIG-2F5-5
Species Reactivity	mouse
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.