

## Purified anti-mouse GM-CSF

Catalog No: #FC1655

Package Size: #FC1655-1 100ug #FC1655-2 1mg

## Description

|                    |                            |
|--------------------|----------------------------|
| Product Name       | Purified anti-mouse GM-CSF |
| Host Species       | rat                        |
| Clonality          | monoclonal                 |
| Clone No.          | s-MP1-22E9                 |
| Species Reactivity | mouse                      |
| Conjugates         | Unconjugated               |
| Storage            | Store at 4°C.              |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.