

Purified anti-mouse IFNAR-1

Catalog No: #FC1660

Package Size: #FC1660-1 100ug #FC1660-2 1mg

Description

Product Name	Purified anti-mouse IFNAR-1
Host Species	mouse
Clonality	monoclonal
Clone No.	s-MAR1-5A3
Species Reactivity	mouse
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.