

Purified anti-mouse IL-1 α

Catalog No: #FC1663

Package Size: #FC1663-1 100ug #FC1663-2 1mg

Description

Product Name	Purified anti-mouse IL-1 α
Host Species	Armenian Hamster
Clonality	monoclonal
Clone No.	s-ALF-161
Species Reactivity	mouse
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.