

Purified anti-human/mouse TGF- β 1, 2, 3

Catalog No: #FC1676

Package Size: #FC1676-1 100ug #FC1676-2 1mg

Description

Product Name	Purified anti-human/mouse TGF- β 1, 2, 3
Host Species	mouse
Clonality	monoclonal
Clone No.	s-1D11.16.8
Species Reactivity	human/mouse
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.