

Purified Syrian Hamster IgG isotype Control

Catalog No: #FC1685

Package Size: #FC1685-1 100ug #FC1685-2 1mg



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description	
Product Name	Purified Syrian Hamster IgG isotype Control
Species Reactivity	Syrian Hamster
Conjugates	Unconjugated
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.