

Histone H3 (L-Lactyl-Lys18) Rabbit mAb

Catalog No: #HL014

Package Size: #HL014-1 100ul

Description

Product Name	Histone H3 (L-Lactyl-Lys18) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB;IP
Species Reactivity	Human;Mouse
Conjugates	Unconjugated
Target Name	Histone H3
Other Names	H3/a, H3/c, H3/d, H3/f, H3/h
Accession No.	P68431
Calculated MW	15 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20Λ C

Application Details

WB 1:500-1:2000; IP 1:20-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.