

Histone H2B (L-Lactyl-Lys15) Rabbit mAb

Catalog No: #HL015

Package Size: #HL015-1 100ul

Description

Product Name	Histone H2B (L-Lactyl-Lys15) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB;IP
Species Reactivity	Human;Mouse
Conjugates	Unconjugated
Target Name	Histone H2B
Other Names	HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF; H2BFG; HIST1H2BG; H2BFA; HIST1H2BI; H2BFK; Histone H2B type 1-C/E/F/G/I; Histone H2B.1 A; Histone H2B.a; H2B/a; Histone H2B.g; H2B/g; Histone H2B.h; H2B/h; Histone H2B.k; H2B/k; Histone H2B.l;
Accession No.	P62807
Calculated MW	14 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20 $\wedge$ C

Application Details

WB 1:500-1:2000; IP 1:20-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.