

Histone H3(Tri-methyl-K9) Mouse Monoclonal Antibody

Catalog No: #HW051

Package Size: #HW051 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Histone H3(Tri-methyl-K9) Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	The Histone H3(Tri-methyl-K9) antibody can detects endogenous Histone H3(Tri-methyl-K9) protein.
Immunogen Description	A synthetic Tri-methylated peptide corresponding to residues surrounding K9 of human histone H3.
Target Name	Histone H3(Tri-methyl-K9)
Modification	Methyl
Accession No.	Swiss-Prot#:P68431
SDS-PAGE MW	15 kDa
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C


Application Details

WB 1:1000-3000


Images



Western blot analysis of Hela cells, using #HW051 diluted at
1) 1:2,000 2) 1:5,000



Western blot analysis of 1) Rat Testis Tissue, 2) Raw264.7, using #HW051 diluted at 1:2000.



1B'B'Input: HeLa Cell Lysate
 2B'B'IP product: IP dilute 1:200
 Western blot analysis: primary antibody : #HW051 1:1,000
 Secondary antibody: Goat anti-Mouse IgG(H+L) 1:10,000

Background

Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.