

Histone H3(Di-methyl-K79) Mouse Monoclonal Antibody (362)

Catalog No: #HW053

Package Size: #HW053 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Histone H3(Di-methyl-K79) Mouse Monoclonal Antibody (362)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	362
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	The Histone H3(Di-methyl-K79) antibody can detects endogenous Histone H3(Di-methyl-K79) protein.
Immunogen Description	A synthetic Di-methylated peptide corresponding to residues surrounding K79 of human histone H3.
Target Name	Histone H3(Di-methyl-K79)
Modification	Methyl
Accession No.	Swiss-Prot#:P68431
SDS-PAGE MW	17kd
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000~1:3000

Images



Western blot analysis of 1) Hela, 2) Raw264.7, 3) Rat Brain Tissue, using #HW053 diluted at 1:2,000.

Background

Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome.

Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.