

NFkB-p65 (Acetyl-Lys310) Polyclonal Antibody

Catalog No: #HW110

Package Size: #HW110-1 50ul #HW110-2 100ul

Description

Product Name	NFkB-p65 (Acetyl-Lys310) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	Acetyl-NFkB-p65 (K310) Polyclonal Antibody detects endogenous levels of NFkB-p65 protein only when acetylated at K310.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human NFkB-p65 around the acetylation site of K310.
Conjugates	Unconjugated
Target Name	NFkB-p65
Modification	Acetyl
Other Names	RELA; NFkB3; Transcription factor p65; Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3
Accession No.	Swiss-Prot: Q04206NCBI Gene ID: 5970
Target Species	human
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000


ELISA: 1/40000

Not yet tested in other applications

Images



Western Blot analysis of HeLa cells using Acetyl-NFκB-p65 (K310) Polyclonal Antibody



Western Blot analysis of customer's sample cells using Acetyl-NFκB-p65 (K310) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.