

HSP 90 (Acetyl-Lys435) Polyclonal Antibody

Catalog No: #HW119

Package Size: #HW119-1 50ul #HW119-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HSP 90 (Acetyl-Lys435) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Acetyl-HSP 90 (K435) Polyclonal Antibody detects endogenous levels of HSP 90 protein only when acetylated at K435.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the human HSP 90 around the acetylation site of K435.
Conjugates	Unconjugated
Target Name	HSP 90
Modification	Acetyl
Other Names	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38
Accession No.	Swiss-Prot: P07900NCBI Gene ID: 3320
Target Species	human
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000


ELISA: 1/20000

Not yet tested in other applications

Images



Western Blot analysis of A549 cells using Acetyl-HSP 90 (K435) Polyclonal Antibody



Western Blot analysis of A549 cells using Acetyl-HSP 90 (K435) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.