

Tubulin α (Acetyl-Lys112) Polyclonal Antibody

Catalog No: #HW131

Package Size: #HW131-1 50ul #HW131-2 100ul

Description

Product Name	Tubulin α (Acetyl-Lys112) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Acetyl-Tubulin α (K112) Polyclonal Antibody detects endogenous levels of Tubulin α protein only when acetylated at K112.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Tubulin α around the acetylation site of K112.
Conjugates	Unconjugated
Target Name	Tubulin α
Modification	Acetyl
Other Names	TUBA1A; TUBA3; Tubulin alpha-1A chain; Alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; Alpha-tubulin ubiquitous
Accession No.	Swiss-Prot: Q71U36/P68363/Q9BQE3NCBI Gene ID: 7846/10376/84790
Target Species	human
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000

ELISA: 1/20000

Not yet tested in other applications

Images



Western Blot analysis of HepG2 cells using Acetyl-Tubulin α (K112) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.