

Histone H3 (Di-Methyl-Lys36) Monoclonal Antibody

Catalog No: #HW167

Package Size: #HW167-1 50ul #HW167-2 100ul


Description

Product Name	Histone H3 (Di-Methyl-Lys36) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Histone H3 (di methyl K36) protein.
Immunogen Type	peptide
Immunogen Description	di methyl peptide at K36 of Histone H3
Conjugates	Unconjugated
Target Name	Histone H3
Modification	Methyl
Other Names	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/
Accession No.	Swiss-Prot: P68431/Q71DI3/P84243NCBI Gene ID: 8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
SDS-PAGE MW	15kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of 1) HeLa, 2) Raw264.7, 3) Rat Brain Tissue, diluted at 1:3,000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.