

Histone H2B (Acetyl-Lys24/25) Polyclonal Antibody

Catalog No: #HW177

Package Size: #HW177-1 50ul #HW177-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Histone H2B (Acetyl-Lys24/25) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Histone H2B when Acetyl occurs at Lys24 or 25)
Immunogen Type	protein
Immunogen Description	Synthetic Acetyl peptide from human protein at AA range: 24/25
Conjugates	Unconjugated
Target Name	Histone H2B (Acetyl-Lys24/25)
Accession No.	Swiss-Prot:Q96A08/P33778/P62807NCBI Gene ID:255626
Calculated MW	15 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western blot analysis of K562 mouse-lung lysate, antibody was diluted at 2000. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.