

Rb (Acetyl-K873/K874) Polyclonal Antibody

Catalog No: #HW189

Package Size: #HW189-1 50ul #HW189-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Rb (Acetyl-K873/K874) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	This antibody detects endogenous levels of Rb (Acetyl-K873/K874). It doesn't react with total protein.
Immunogen Type	protein
Immunogen Description	Synthesized Acetyl peptide derived from human Rb. at AA range: K873/872
Conjugates	Unconjugated
Target Name	Rb (Acetyl-K873/K874)
Other Names	Retinoblastoma-associated protein (p105-Rb) (pRb) (Rb) (pp110)
Accession No.	Swiss-Prot:P06400NCBI Gene ID:5925
Calculated MW	106 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,mouse-heart 3,mouse-brain cells using primary antibody diluted at 1:1000(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B1°C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.