

Histone H4 (Mono Methyl Lys79) Polyclonal Antibody

Catalog No: #HW191

Package Size: #HW191-1 50ul #HW191-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Histone H4 (Mono Methyl Lys79) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Histone H4 (Mono Methyl Lys79) protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Histone H4 (Mono Methyl Lys79)
Conjugates	Unconjugated
Target Name	Histone H4
Other Names	HIST1H4A; H4/A; H4FA; HIST1H4B; H4/I; H4FI; HIST1H4C; H4/G; H4FG; HIST1H4D; H4/B; H4FB; HIST1H4E; H4/J; H4FJ; HIST1H4F; H4/C; H4FC; HIST1H4H; H4/H; H4FH; HIST1H4I; H4/M; H4FM; HIST1H4J; H4/E; H4FE; HIST1H4K; H4/D; H4FD; HIST1H4L; H4/K; H4FK; HIST2H4A; H4/N; H4F2; H4FN; HIST2H4; HIST2H4B; H4/O; H4FO;
Accession No.	Swiss-Prot:P62805Gene ID:121504
Calculated MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:1000

Images



Western blot analysis of 1) Hela, 2) 3T3, diluted at 1:2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.