

Histone H1 (Mono Methyl Lys25) Polyclonal Antibody

Catalog No: #HW192

Package Size: #HW192-1 50ul #HW192-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Histone H1 (Mono Methyl Lys25) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Histone H1 (Mono Methyl Lys25) protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Histone H1 (Mono Methyl Lys25)
Conjugates	Unconjugated
Target Name	Histone H1
Other Names	HIST1H1A; H1F1; Histone H1.1; Histone H1a; HIST1H1C; H1F2; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1; HIST1H1D; H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Histone H1.4; Histone H1b; Histone H1s-4
Accession No.	Swiss-Prot:Q02539/P16403/P16402/P10412Gene ID:3024/3006/3007/3008
Calculated MW	17-25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:1000

Images



Western blot analysis of 1) Hela, 2) 3T3, diluted at 1:2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.