

Histone H2B (Mono Methyl Lys5) Polyclonal Antibody

Catalog No: #HW193

Package Size: #HW193-1 50ul #HW193-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Histone H2B (Mono Methyl Lys5) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Histone H2B (Mono Methyl Lys5) protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Histone H2B (Mono Methyl Lys5)
Conjugates	Unconjugated
Target Name	Histone H2B
Other Names	HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF; H2BFG; HIST1H2BG; H2BFA; HIST1H2BI; H2BFK; Histone H2B type 1-C/E/F/G/I; Histone H2B.1
Accession No.	Swiss-Prot:Q96A08/P33778/P62807Gene ID:255626/3018/3017
Calculated MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:1000

Images



Western blot analysis of 1) HeLa, 2) 3T3, diluted at 1:2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.