

Histone H4 (Mono Methyl Arg3) Polyclonal Antibody

Catalog No: #HW194

Package Size: #HW194-1 50ul #HW194-2 100ul


Description

Product Name	Histone H4 (Mono Methyl Arg3) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG1
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous Histone H3 (Mono Methyl Arg3) protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Histone H4 (Mono Methyl Arg3)
Conjugates	Unconjugated
Target Name	Histone H4
Other Names	HIST1H4A; H4/A; H4FA; HIST1H4B; H4/I; H4FI; HIST1H4C; H4/G; H4FG; HIST1H4D; H4/B; H4FB; HIST1H4E; H4/J; H4FJ; HIST1H4F; H4/C; H4FC; HIST1H4H; H4/H; H4FH; HIST1H4I; H4/M; H4FM; HIST1H4J; H4/E; H4FE; HIST1H4K; H4/D; H4FD; HIST1H4L; H4/K; H4FK; HIST2H4A; H4/N; H4F2; H4FN; HIST2H4; HIST2H4B; H4/O; H4FO;
Accession No.	Swiss-Prot:P62805Gene ID:121504
Calculated MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-1:2000

Images



Western blot analysis of extracts from Hela cells, 1:2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.