

α -tubulin (Acetyl Lys40) Monoclonal Antibody(5H5)

Catalog No: #HW204

Package Size: #HW204-1 50ul #HW204-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	α -tubulin (Acetyl Lys40) Monoclonal Antibody(5H5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5H5
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	α -tubulin (Acetyl Lys40) Mouse Monoclonal Antibody(5H5) detects endogenous levels of α -tubulin (Acetyl Lys40)
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	α -tubulin
Accession No.	Swiss-Prot:P68363Gene ID:10376
Calculated MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:100-1:200

Images



Western blot analysis of extracts from HeLa cells, untreated (-) or treated with TSA (1 μ M, 18 hr; +), using Acetyl- α -tubulin(Lys40) Mouse mAb 1:2000.



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using α -tubulin(Acetyl Lys40) Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.