

Acetyl NF kB P65B£B"K314/K315) Mouse Monoclonal Antibody(5G11)

Catalog No: #HW207

Package Size: #HW207-1 50ul #HW207-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Acetyl NF kB P65B£B"K314/K315) Mouse Monoclonal Antibody(5G11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5G11
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Acetyl NF kB P65B£B"K314/K315) protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Acetyl NF kB P65B£B"K314/K315)
Conjugates	Unconjugated
Target Name	NF kB P65
Other Names	Transcription factor p65 (Nuclear factor NF-kappa-B p65 subunit) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3)
Accession No.	Swiss-Prot:Q04206Gene ID:5970
Calculated MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl NF kB P65(K314/K315) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Acetyl NF κB P65(K314/K315) Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.