

Rabbit IgG Isotype Control PE-Cy5 Conjugated

Catalog No: #IC001PE-Cy5

Package Size: #IC001PE-Cy5-1 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit IgG Isotype Control PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB, ELISA, FCM, IHC-P, IHC-F, IF(IHC-P), IF(IHC-F), IF(ICC)
Conjugates	PE-Cy5
Concentration	1ug/ul
Storage	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide. Store at -20°C for 12 months.

Application Details

WB Dilution:1:300-1000
ELISA Dilution:1:500-1000
FCM Dilution:1:20-100
IHC-P Dilution:1:200-400
IHC-F Dilution:1:100-500
IF(IHC-P) Dilution:1:50-200
IF(IHC-F) Dilution:1:50-200
IF(ICC) Dilution:1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.