

FITC mouse IgG1,k Isotype Control

Catalog No: #IC002FITC

Package Size: #IC002FITC 50tests

Description

Product Name	FITC mouse IgG1,k Isotype Control
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1,k
Purification	Purified by Protein A.
Applications	FC IF
Conjugates	FITC
Concentration	10 $\times$ OI/Test
Storage	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA,Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Background

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.