

Goat Anti-Rabbit IgG Secondary Antibody AF488
Conjugated

Catalog No: #L3016

Package Size: #L3016-1 1.0 ml



Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody AF488 Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	AF488
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:Whole IgG

most applications: 1:100 - 1:800

Published Papers

el at., Gonadal development and molecular analysis revealed the critical window for sex differentiation, and E2 reversibility of XY-male spotted scat, Scatophagus argus., , (2021)
[PMID:](#)

el at., A Novel Airway-Organoid Model Based on a Nano-Self-Assembling Peptide: Construction and Application in Adenovirus Infection StudiesIn Int J NanomedicineOn2023 Sep 15byYun-E Xu , Di-Shu Ao et al..PMID:37727653, , (2023)
[PMID:37727653](#)

el at., PTBP1 knockdown impairs autophagy flux and inhibits gastric cancer progression through TXNIP-mediated oxidative stress. In Cell Mol Biol Lett on 2024 Aug 17 by Shimin Wang, Xiaolin Wang,et al..PMID:39153986, , (2024)
[PMID:39153986](#)

Weiping Xu;Shujuan Yi;Jie Liu;Yuyan Jiang;Jiaguo Huang el at., Nitrile-aminothiol bioorthogonal near-infrared fluorogenic probes for ultrasensitive in vivo imaging., , (2025)
[PMID:39747031](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.