

Goat Anti-Human IgG Secondary Antibody HRP Conjugated

Catalog No: #L3022

Package Size: #L3022-1 1ml



Description

Product Name	Goat Anti-Human IgG Secondary Antibody HRP Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	IgG (H+L)
Conjugates	Horseradish Peroxidase
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:Whole IgG
IHC/ICC: 1:500 - 1:5,000
ELISA: 1:5,000 - 1:50,000
Western blot: 1:5,000 - 1:50,000

Published Papers

Lili Zhao;Jiahao Du;Yuxin Jin;Ying Hu;Suqin Zhang;Biao Yang;Chenglin Sun;Yunxiao Zhao;Xinxin Ge;Rong Yan;Chunliang Liu;Renping Hu;Kesheng Dai et al., A novel antibody against GPIIb α inhibits platelet function and thrombosis without increasing bleeding., , (2025)
[PMID:39924011](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.