

Goat Anti-Mouse IgG Secondary Antibody Rhodamine (RRX) Conjugated

Catalog No: #L30315

Package Size: #L30315-1 1.0 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody Rhodamine (RRX) Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	Rhodamine Red-X (RRX)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.