

# Goat Anti-Mouse IgG Secondary Antibody PE Conjugated

Catalog No: #L30317

Package Size: #L30317-1 1.0 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                               |
|--------------------|-----------------------------------------------------------------------------------------------|
| Product Name       | Goat Anti-Mouse IgG Secondary Antibody PE Conjugated                                          |
| Host Species       | Goat                                                                                          |
| Clonality          | Polyclonal                                                                                    |
| Species Reactivity | Ms                                                                                            |
| Specificity        | IgG (H+L)                                                                                     |
| Conjugates         | R-Phycoerythrin                                                                               |
| Formulation        | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, |
| Storage            | Store at -20 °C/1 year                                                                        |

## Application Details

|                                  |
|----------------------------------|
| Antibody Format:F(ab') Fragment  |
| most applications: 1:100 - 1:200 |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.