

Goat Anti-Pig IgG (H&L) HRP

Catalog No: #L35008

Package Size: #L35008-1 100 ul #L35008-2 500 ul #L35008-3 1 ml

Description

Product Name	Goat Anti-Pig IgG (H&L) HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	P
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Pig IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Pig IgG
Conjugates	HRP
Concentration	1 mg/ml
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.