

# Rabbit Anti-Horse IgG (H&L) HRP

Catalog No: #L35019

Package Size: #L35019 100 ul

## Description

|                       |                                                                                                                                                          |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Rabbit Anti-Horse IgG (H&L) HRP                                                                                                                          |
| Host Species          | Rabbit                                                                                                                                                   |
| Clonality             | Polyclonal                                                                                                                                               |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.                                         |
| Applications          | E, WB, IH, IC                                                                                                                                            |
| Species Reactivity    | Hs                                                                                                                                                       |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Horse IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Horse IgG                                                                                                                                                |
| Conjugates            | HRP                                                                                                                                                      |
| Concentration         | 1mg/ml                                                                                                                                                   |
| Formulation           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin                                                        |
| Storage               | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                         |

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.