

## Mouse Anti-Chicken IgY (H&amp;L) HRP

Catalog No: #L35029

Package Size: #L35029-1 100 ul #L35029-2 500 ul #L35029-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                                            |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Mouse Anti-Chicken IgY (H&L) HRP                                                                                                                           |
| Host Species          | Mouse                                                                                                                                                      |
| Clonality             | Monoclonal                                                                                                                                                 |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.                                           |
| Applications          | E, WB, IH, IC                                                                                                                                              |
| Species Reactivity    | C                                                                                                                                                          |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Chicken IgY                                                                                                                                                |
| Conjugates            | HRP                                                                                                                                                        |
| Formulation           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide                                                     |

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.