

Rabbit Anti-Human IgG+IgM+IgA Biotin

Catalog No: #L35054

Package Size: #L35054-1 100 ul #L35054-2 500 ul #L35054-3 1 ml

Description

Product Name	Rabbit Anti-Human IgG+IgM+IgA Biotin
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	H
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG IgM IgA. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human IgG+IgM+IgA
Conjugates	Biotin
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.