

Rabbit Anti-Human IgG (H&L) FITC

Catalog No: #L35084

Package Size: #L35084-1 100 ul #L35084-2 500 ul #L35084-3 1 ml



**SAB**<sup>®</sup>  
Signalway Antibody

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                       |                                                                                                                                                          |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Rabbit Anti-Human IgG (H&L) FITC                                                                                                                         |
| Host Species          | Rabbit                                                                                                                                                   |
| Clonality             | Polyclonal                                                                                                                                               |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.                                         |
| Applications          | IF, FC                                                                                                                                                   |
| Species Reactivity    | H                                                                                                                                                        |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgG                                                                                                                                                |
| Conjugates            | FITC                                                                                                                                                     |
| Formulation           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide                                                   |

Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.