

Rabbit Anti-Human IgG+IgM+IgA FITC

Catalog No: #L35086

Package Size: #L35086-1 100 ul #L35086-2 500 ul #L35086-3 1 ml



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Human IgG+IgM+IgA FITC
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	H
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG IgM IgA. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human IgG+IgM+IgA
Conjugates	FITC
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.