

## Goat Anti-Human IgG+IgM+IgA - AF405

Catalog No: #L35145

Package Size: #L35145-1 100 ul #L35145-2 500 ul #L35145-3 1 ml

## Description

|                       |                                                                                                                                                                  |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Goat Anti-Human IgG+IgM+IgA - AF405                                                                                                                              |
| Host Species          | Goat                                                                                                                                                             |
| Clonality             | Polyclonal                                                                                                                                                       |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.                                                |
| Applications          | E, IH, IF, FC                                                                                                                                                    |
| Species Reactivity    | H                                                                                                                                                                |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG IgM IgA. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgG+IgM+IgA                                                                                                                                                |
| Conjugates            | AF Dye                                                                                                                                                           |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide                                                  |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.