

## Goat Anti-Mouse IgM mu chain - AF568

Catalog No: #L35320

Package Size: #L35320-1 100 ul #L35320-2 500 ul #L35320-3 1 ml

## Description

|                       |                                                                                                                                                                   |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Goat Anti-Mouse IgM mu chain - AF568                                                                                                                              |
| Host Species          | Goat                                                                                                                                                              |
| Clonality             | Polyclonal                                                                                                                                                        |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.                                                 |
| Applications          | E, IH, IF, FC                                                                                                                                                     |
| Species Reactivity    | M                                                                                                                                                                 |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgM mu chain. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Mouse IgM mu chain                                                                                                                                                |
| Conjugates            | AF Dye                                                                                                                                                            |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide                                                   |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.