

# Goat Anti-Human IgG Fc - AF594

Catalog No: #L35358

Package Size: #L35358-1 100 ul #L35358-2 500 ul #L35358-3 1 ml

## Description

|                       |                                                                                                                                                             |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Goat Anti-Human IgG Fc - AF594                                                                                                                              |
| Host Species          | Goat                                                                                                                                                        |
| Clonality             | Polyclonal                                                                                                                                                  |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.                                           |
| Applications          | E, IH, IF, FC                                                                                                                                               |
| Species Reactivity    | H                                                                                                                                                           |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG Fc. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgG Fc                                                                                                                                                |
| Conjugates            | AF Dye                                                                                                                                                      |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide                                             |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.