

Rabbit Anti-Sheep IgG (H&L) - AF594

Catalog No: #L35384

Package Size: #L35384-1 100 ul #L35384-2 500 ul #L35384-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Sheep IgG (H&L) - AF594
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	S
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Sheep IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Sheep IgG
Conjugates	AF Dye
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.