

## Rabbit Anti-Goat IgG (H&amp;L) - AF633

Catalog No: #L35415

Package Size: #L35415-1 100 ul #L35415-2 500 ul #L35415-3 1 ml

## Description

|                       |                                                                                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Rabbit Anti-Goat IgG (H&L) - AF633                                                                                                                      |
| Host Species          | Rabbit                                                                                                                                                  |
| Clonality             | Polyclonal                                                                                                                                              |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.                                       |
| Applications          | E, IH, IF, FC                                                                                                                                           |
| Species Reactivity    | G                                                                                                                                                       |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Goat IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Goat IgG                                                                                                                                                |
| Conjugates            | AF Dye                                                                                                                                                  |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide                                         |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.