

Goat Anti-Chicken IgY (H&L) - AF647

Catalog No: #L35441

Package Size: #L35441 100 ul

Description

Product Name	Goat Anti-Chicken IgY (H&L) - AF647
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	C
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Chicken IgY
Conjugates	AF Dye
Concentration	1 mg/ml
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.